

Fractional Laser Machine - Official Clinical Overview & Datasheet

EXECUTIVE SUMMARY

This document presents the clinical and technical overview of the next-generation Fractional Laser System, a non-ablative and ablative hybrid platform engineered for dermatological resurfacing, scar revision, and photoaging correction. The device is clinically positioned for high-volume medical spas, dermatology clinics, and plastic surgery centers seeking a versatile solution that bridges the gap between aggressive CO2 lasers and low-power non-ablative devices. With a proprietary pulse stacking technology and real-time skin impedance monitoring, the system delivers controlled micro-injury zones (MTZs) while minimizing thermal spread, thereby reducing patient downtime from 14 days to an average of 3-5 days.

The primary value proposition centers on three pillars: superior clinical efficacy for deep rhytides and atrophic scars, an integrated epidermal protection system allowing treatment on Fitzpatrick skin types IV-VI, and a high-ROI operational model (consumable-free handpieces, 50,000+ pulse durability, and 30-second cartridge swap). For clinic owners, the device achieves an average of 4.2 treatments per complete scar revision protocol, yielding cost-per-procedure reductions of 40% compared to fractional CO2 alternatives.



CLINICAL ARCHITECTURE & DESIGN

The hardware architecture employs a dual-source resonator cavity using air-cooled, hermetically sealed fractional laser bars imported from Germany. Unlike conventional systems that rely on single flashlamps, this design incorporates two independent 300W diode stacks (pump sources) driving a Nd:YAG or Er:YAG crystal—selectable by the clinician via the handpiece interface. Each laser bar features gold-plated heat sinks and real-time temperature telemetry with $\pm 0.5^{\circ}$ C accuracy, ensuring wavelength stability across continuous 8-hour operation. The fluid dynamics subsystem includes a magnetic-drive water pump (MTBF 10,000 hours), dual 120mm high-static pressure fans, and a 2.5L corrosion-resistant reservoir with automated deaeration, eliminating manual priming.

Advanced epidermal cooling is achieved via a three-stage cascade: (1) contact

sapphire tip pre-cooling to -2°C within 0.2 seconds of skin contact, (2) dynamic cold air pulse synchronized to each fractional pulse, and (3) post-treatment cryogen spray with adjustable duration (10-80ms). The handpiece itself contains a MEMS accelerometer that disables laser emission if motion exceeds 20mm/s, preventing overlap burns. The smart UI includes a 15.6-inch capacitive touchscreen with glove-friendly interface, preset protocols for 12 indications (e.g., melasma, surgical scars, perioral wrinkles), and a live thermal camera overlay showing epidermal temperature gradients.

KEY INDICATIONS & CAPABILITIES

- Fractional Resurfacing (Deep & Superficial): Selectable depth control from $50\mu\text{m}$ (superficial) to 1.5mm (deep dermis) via adjustable pulse energy (5-120mJ/microbeam). Achieves 40-60% collagen contraction per session with 15-20% treatment density, validated for Fitzpatrick II-V.
- Acne & Trauma Scar Revision: Proprietary "stack and pause" algorithm delivers 3-5 stacked pulses at 15Hz into each micro-channel, breaking fibrotic septae while preserving intervening tissue. Clinical study (n=87) showed 68% improvement in Goodman and Baron scar scores after 3 sessions.
- Pigmentation & Melasma Clearance: Low-fluence fractional mode (5-10mJ, 10% density) combined with dual-wavelength option (1540nm non-ablative) bypasses the basal layer to target dermal melanophages. Post-inflammatory hyperpigmentation rate 8% in Fitzpatrick IV (traditional CO₂: 32%).

- Smart UI & Safety Interlocks: Multi-factor authentication (RFID key + password) for user level management; real-time contact quality assessment with resistance monitoring (20-200 kOhm acceptable range); automatic fluence derating by 15% per 1mm of skin curvature detected.
- Versatile Handpiece Ecosystem: Snap-on disposable tips for hygiene (20mm, 12mm, 8mm spot sizes) plus reusable autoclavable handpiece body. Tip ID chip auto-configures energy limits and pulse patterns (square, hexagonal, or random microbeam arrays).

COMPLIANCE & STANDARDS

This fractional laser system is manufactured under ISO 13485:2021 certified production and quality management system. It holds Medical CE certification (Class IIb) under EU MDR 2017/745, confirmed by notified body verification of essential requirements annex I. For international markets, the device complies with IEC 60601-1 (medical electrical equipment safety), IEC 60825-1 (laser product safety, Class 4), and IEC 60601-2-22 (particular requirements for surgical, cosmetic, therapeutic laser equipment). FDA 510(k) clearance has been granted for the indication "treatment of wrinkles, pigmented lesions, and surgical scars" (Kxxxxxx). Additional certifications include RoHS 3 compliance, REACH SVHC declaration, and CB Scheme certification for 50+ countries.

TECHNICAL SPECIFICATIONS

The following parameters define the baseline performance envelope of the Fractional Laser System under standard operating conditions (23°C ambient, 40% relative humidity). All values are factory-calibrated and field-verifiable using the included power meter handpiece.

| Parameter | Specification |
|----------------------------|---|
| Laser Type / Wavelength | Ablative: 2940nm Er:YAG / Non-ablative: 1540nm or 1550nm (selectable via handpiece) |
| Fractional Beam Spot Size | 100-300µm (individual microbeam), 10-20mm (overall treatment area) |
| Pulse Energy per Microbeam | 5-120mJ (adjustable in 5mJ increments) |
| Maximum Treatment Density | 5% to 40% (selectable microbeam spacing patterns) |
| Spot Size (Handpiece tip) | 8mm, 12mm, 20mm (interchangeable, RFID-coded) |
| Cooling System | Contact sapphire (-2°C) + Cryogen spray (10-80ms) + Forced air (0-30 L/min) |

| | |
|------------------------|--|
| Repetition Rate | 1-30 Hz (continuous or burst mode) |
| Aiming Beam | 635nm red diode, power <1mW (Class 2M) |
| Display & UI | 15.6" capacitive touch, 1920x1080, real-time skin resistance + thermal overlay |
| Dimensions (Main Unit) | 450mm (W) x 500mm (D) x 1100mm (H) |
| Weight | 42 kg (excluding articulating arm) |
| Power Supply | 100-240V AC, 50/60Hz, 1500VA typical |

