

Q-switched Nd:YAG Laser - Official Clinical Overview & Datasheet

EXECUTIVE SUMMARY

The Q-switched Nd:YAG laser represents the gold-standard platform for professional pigment clearance, tattoo removal, and multi-application cutaneous therapy. Designed for high-volume dermatology clinics and medical spas, this system delivers nanosecond-domain pulse energies that selectively photoacoustically disrupt dermal and epidermal chromophores with minimal thermal diffusion. Its dual-wavelength architecture (1064nm & 532nm) enables precise targeting of deep and superficial pigments, ranging from professional tattoos to melasma, nevus of Ota, and solar lentigines.

Clinically positioned as a high-ROI capital investment, the platform combines therapeutic versatility with low consumable costs and extended component lifespan. Primary value propositions include: pain-manageable treatments without general anesthesia, superior clearance rates over alternative modalities, and built-in epidermal protection via a high-speed cryogen or contact cooling interface. For practice owners, the device reduces per-procedure time, expands service menu offerings, and ensures regulatory compliance for medical reimbursable indications.



CLINICAL ARCHITECTURE & DESIGN

The system integrates a military-grade, air- and water-cooled solid-state Nd:YAG resonator with a premium Q-switch module (either Pockels cell or passive saturable absorber), delivering pulse durations in the 5–20 nanosecond range. The laser cavity incorporates imported optically pumped semiconductor or flashlamp excitation units, ensuring $\pm 3\%$ pulse-to-pulse energy stability over an eight-hour clinical shift. A dual-light guide articulated arm (or high-transmission fiber delivery) terminates in an ergonomic, autoclaveable handpiece with interchangeable spot sizes (1.5mm to 8mm round or square-top profiles).

Advanced epidermal cooling is achieved via a configurable closed-loop chiller with a sapphire contact window (0°C to +4°C surface temperature) synchronized

to the Q-switched emission. An optional R404a cryogen spray unit (40–80ms pulses, adjustable pre/post delay) provides second skin layer protection for high-fluence treatments. The power supply uses redundant insulated-gate bipolar transistors (IGBTs) and a self-diagnosing water pump with flow sensor and particulate filtration to prevent overheating and guarantee 24/7 operation under duty cycles up to 3 pulses per second.

KEY INDICATIONS & CAPABILITIES

- Dual-wavelength photoacoustic targeting: 1064nm (deep: black/dark blue/green tattoos, dermal melanin, leg veins) and 532nm (superficial: red/orange/brown tattoos, epidermal pigmentation, freckles) with automatic wavelength switching in the handpiece.
- Adjustable pulse energy and spot geometry: energy density (fluence) from 0.5 to 12 J/cm² (532nm) and 1 to 16 J/cm² (1064nm); six round spot sizes plus a 2x8mm fractional-like rectangular mode for textural improvement.
- Smart clinical user interface: 10.4-inch capacitive touchscreen preloaded with over 80 preset protocols (Fitzpatrick I–VI, tattoo colors, lesion types) that dynamically set fluence, repetition rate (1–10 Hz), and cooling parameters; real-time skin impedance monitoring triggers automatic fluence reduction.
- Large-mode applicator for high-speed treatments: 6mm and 8mm handpiece tips enable removal of large benign pigmented lesions (e.g., nevus spilus, café-au-lait) at rates up to 2 cm² per second, reducing procedure time by up to 40%

versus traditional 2–3mm spots.

- Integrated laser safety and dosimetry: built-in energy meter with closed-loop feedback, audible countdown, and a dual-footswitch (standby/fire) preventing accidental emission; patient and operator ophthalmic shields included.

COMPLIANCE & STANDARDS

The Q-switched Nd:YAG laser system carries Medical CE Mark (Class 4 active device) under Directive 2017/745/EU, FDA 510(k) clearance for soft tissue incision, tattoo removal, benign pigmented lesion ablation, and wrinkle treatment, ISO 13485:2016 (production quality system), and IEC 60825-1:2014 (laser product safety class 4). Each unit undergoes factory power output validation, MPE (maximum permissible exposure) testing, and includes a technical file for UKCA, NMPA (China), KFDA, and ANVISA registrations. A post-market surveillance plan and full risk management file (ISO 14971) are available to buyers.

TECHNICAL SPECIFICATIONS

The following parameters represent nominal performance values measured at the handpiece tip under calibrated clinical conditions. All values are subject to standard manufacturing tolerance ($\pm 5\%$) and are certified at system delivery.

Parameter	Specification
Laser Type / Wavelengths	Q-switched Nd:YAG / 1064 nm & 532 nm (frequency-doubled)
Pulse Duration (FWHM)	5 – 10 ns (1064 nm) / ≤ 8 ns (532 nm)
Peak Power	> 20 MW (per pulse at 1064 nm, 1000 mJ)
Spot Sizes	1.5, 2, 3, 4, 5, 6, 8 mm round; 2x8 mm fractional rectangular
Fluence Range (1064 nm)	1 – 16 J/cm ² (dependent on spot size)
Fluence Range (532 nm)	0.5 – 12 J/cm ² (dependent on spot size)
Repetition Rate	Single shot – 10 Hz (adjustable in 0.5 Hz steps)
Cooling System	Contact sapphire (0 °C to +4 °C) + optional R404a cryogen spray
Display & Interface	10.4 inch capacitive touchscreen, preset protocols, energy feedback
Electrical / Dimensions	110-240 VAC, 50/60 Hz, 1000W; 45 kg (99 lb); 450x500x950 mm (WxDxH)

