

Deep Tissue RF Therapy Workstation - Clinical Architecture & Performance
Reference Manual

DEEP TISSUE RF THERAPY WORKSTATION: CLINICAL ARCHITECTURE &
PERFORMANCE REFERENCE MANUAL

EXECUTIVE SUMMARY

The Deep Tissue RF Therapy Workstation represents a paradigm shift in non-invasive volumetric heating technology. Unlike legacy monopolar or bipolar systems that deliver superficial thermal gradients, this platform utilizes a proprietary Synchronized Return-Path Electrode Array (SRPEA) to achieve uniform dermal and subdermal heating up to 20mm depth. The system is engineered for dual-clinic utility: Class I medical device for tissue tightening and Class IIa for temporary pain relief in musculoskeletal indications. This reference document details the hardware topology, thermal safety envelope, and verified treatment parameters for procurement and clinical engineering review.



CLINICAL ARCHITECTURE & DESIGN

The Workstation comprises three core subsystems: a 2.0 MHz Class-E resonant solid-state generator, a closed-loop peristaltic thermoelectric cooling module, and a multi-articulating handpiece with real-time contact impedance monitoring. The RF energy is delivered via capacitive coupling through a 2.5 cm² treatment window, with the return path distributed across two independent grounding pads to eliminate hot spots. A proprietary Adaptive Thermal Feedback loop samples skin surface temperature 1000 times per second via an infrared microbolometer, modulating power output to maintain an exact setpoint of 39°C to 44°C for neocollagenesis without epidermal injury.

KEY INDICATIONS & CAPABILITIES

Clinically validated for Fitzpatrick skin types I-VI. Primary aesthetic indications: reduction of mild-to-moderate facial wrinkles (Fitzpatrick Wrinkle Scale 2-4), circumferential abdominal laxity, and post-partum dermal remodeling. Secondary medical indications: adjunctive treatment of chronic nonspecific low back pain and knee osteoarthritis pain. The workstation delivers a maximum RF power of 120W continuous wave, pulse-mode fractional heating at 50% duty cycle (10Hz), and a unique Thermistor-Controlled Ramping protocol for anesthetic-free treatments. Average session time: 25 minutes (full face) to 40 minutes (abdomen + flanks).

COMPLIANCE & STANDARDS

The device carries CE MDR (Class IIa), UKCA, and South Korea MFDS clearance. FDA 510(k) clearance under product code OLP (Electrosurgical, Cut & Coagulation & Accessories) for tissue coagulation and tightening. Safety certified to IEC 60601-1 (3.1 Ed.), IEC 60601-2-2 (High Frequency Surgical Equipment), and IEC 60601-1-11 (Home Healthcare Environment). Electromagnetic compatibility per CISPR 11 Group 1 Class A. The cooling subsystem meets ISO 80601-2-35 for thermal safety. All patient-facing components are biocompatible per ISO 10993-5 (cytotoxicity) and ISO 10993-10 (irritation).

TECHNICAL SPECIFICATIONS

RF GENERATOR TYPE: Solid-state, Class-E resonant, software-controlled impedance matching

OUTPUT FREQUENCY: 2.0 MHz \pm 5%

MAXIMUM OUTPUT POWER: 120W continuous, 150W peak-pulse (50ms pulse width)

IMPEDANCE RANGE: 50 Ω to 1000 Ω with automatic matching (<50ms lock time)

TREATMENT MODES: Continuous volumetric, Fractionated pulse (1-20Hz), Thermo-ramping (0.5°C/sec)

THERMAL MONITORING: Non-contact IR sensor, 8-14 μ m spectral range, accuracy \pm 0.3°C

SKIN TEMPERATURE LIMITS: 39°C to 44°C (therapeutic), >44°C audio-visual alarm + automatic power cut

COOLING SYSTEM: Sealed recirculating liquid chiller with Peltier pre-cooler; setpoint 5°C to 15°C

HANDPIECE WEIGHT: 280g (active tip included), ergonomic rotary grip

DISPLAY: 10.1" capacitive touchscreen, 1280x800 resolution, anti-glare coating

POWER SUPPLY: Universal 100-240VAC, 50/60Hz, 600VA max

DIMENSIONS (W x D x H): 380mm x 420mm x 980mm (on rolling stand)

WEIGHT: 28kg (device only), 32kg (including stand and accessory tray)

Parameter	Specification
Therapeutic Modality	Deep Tissue Radiofrequency (2.0 MHz synchronous return-path)
Maximum Output Power	120W continuous / 150W pulsed (50ms, 10Hz)
Thermal Safety Envelope	39 ° C - 44 ° C (IR microbolometer feedback, ±0.3°C)
Active Treatment Area	2.5 cm ² (capacitive coupling tip)
Cooling System	Closed-loop liquid chiller + Peltier, 5°C to 15°C setpoint
Duty Cycle / Pulse Modes	Continuous, 10-90% fractionated (1-20Hz), Thermo-ramping
Impedance Adaptation	50Ω to 1000Ω, software auto-match, <50ms lock
Patient Skin Types	Fitzpatrick I - VI (phototypes tested per ISO 24444)
Medical Compliance	CE MDR Class IIa, FDA 510(k) OLP, UKCA, MFDS
Power / Grid Requirements	100-240VAC, 50/60Hz, 600VA, IEC C14 inlet

CLINICAL PROTOCOLS

All treatments require a signed informed consent and a test spot at 38°C for 30 seconds to rule out nickel or thermal urticaria. For facial applications: apply conductive gel, use the gliding technique at 2-3 cm/s with target temperature 42°C \pm 1°C, total fluence per pass 12-15 kJ, 3 passes per session. For body contouring: use stamping technique (3 seconds stationary, 1 second rest) over 10x10cm grids, target temperature 43°C \pm 0.5°C, 8-10 kJ per grid, 2 sessions separated by 21 days. Post-treatment: apply soothing hydrogel; erythema typically resolves in 30-90 minutes. Expected collagen remodeling visible at 8-12 weeks. Maximum three sessions per anatomical area per year.

