

Monopolar RF Skin Tightening Machine - Official Clinical Overview & Technical Datasheet

DEVICE IDENTIFICATION & ROLE IN DERMATOLOGY

The Monopolar RF Skin Tightening Machine represents a Class IIa medical device engineered for non-invasive dermal remodeling. Operating on the principle of volumetric resistive heating, the device delivers controlled deep dermal coagulation while preserving the epidermal layer. This platform is indicated for mild-to-moderate skin laxity reduction on the face, neck, abdomen, and extremities.



INTERNAL HARDWARE TOPOLOGY

The system integrates a high-efficiency 1.0 MHz – 4.0 MHz tunable RF

generator coupled with a closed-loop impedance monitoring circuit. The mainboard architecture includes an isolated power supply unit (PSU), a capacitive touch interface controller, and dual-channel temperature feedback processors. Energy is transmitted through a single-pole active electrode handpiece, with a large dispersive return pad ensuring complete circuit closure. An onboard microcontroller performs 1,000+ impedance readings per second to modulate energy output in real-time, preventing localized hotspots.

EPIDERMAL PROTECTION MECHANISMS

To ensure patient safety across Fitzpatrick skin types I-VI, the device employs a contact-cooled sapphire tip (0 ° C to +4 ° C) integrated into the treatment handpiece. This active cooling system creates a thermal gradient: the epidermis remains below 42 ° C while the target dermal layer reaches therapeutic temperatures of 55°C to 65°C. Additional safety layers include: (1) automatic energy shut-off if skin impedance rises above 180 ohms, (2) real-time temperature monitoring via an infrared thermopile, and (3) a motion detection sensor that halts RF output if the handpiece remains stationary for >0.5 seconds.

TREATMENT ADVANTAGES

- Immediate collagen contraction (primary effect) followed by neocollagenesis over 3-6 months
- Uniform volumetric heating up to 5mm – 8mm penetration depth
- Zero post-treatment downtime with minimal patient discomfort
- Applicable to periorbital, submental, and body contouring zones
- No consumable gels required beyond standard conductive gel

SPECIFICATION MATRIX

Parameter	Specification
RF Frequency	1.0 MHz \pm 10% (Tunable to 4.0 MHz optional)
Maximum Output Power	200W (Continuous wave)
Impedance Range	20 – 180 Ω (Auto-adjusting)
Cooling System	TEC + Sapphire contact cooling (0°C to +20°C adjustable)
Handpiece Active Tip Area	3.0 cm ² (Standard); 1.5 cm ² (Periorbital optional)
Penetration Depth	5-8 mm (At 1.0 MHz)
Temperature Monitoring	Dual IR sensors (\pm 0.5°C accuracy)
Display	8-inch Capacitive Touchscreen, 1280x800 resolution

Power Supply	AC 110-240V, 50/60Hz, 600VA max
Unit Dimensions (WxDxH)	380mm x 420mm x 180mm
Weight	12.5 kg (Main console only)

REGULATORY COMPLIANCE

The device holds CE Mark (MDR 2017/745) and 510(k) clearance for over-the-counter secondary sale in specific markets. Compliance includes IEC 60601-1 (Medical Electrical Equipment Safety), IEC 60601-2-2 (High Frequency Surgical Equipment), and IEC 60601-1-11 (Home Healthcare Environment). The system is manufactured under ISO 13485:2016 certified facilities, with full electromagnetic compatibility (EMC) testing per CISPR 11 Group 1 Class A.



CLINICAL PROTOCOLS (EXCERPT)

- Pre-treatment: Remove all metal jewelry and apply 2-3mm layer of conductive gel.
- Parameter initiation: 15-30W for facial applications; 30-50W for body.
- Technique: slow circular passes (0.5-1 cm/sec) with continuous skin contact.
- Endpoint: mild erythema (1-2/10) and surface temperature of 39°C-41°C.
- Session duration: 20-40 minutes per anatomical area.
- Protocol: 3-6 sessions spaced 4 weeks apart; maintenance every 6 months.