

## SHR Hair Removal Machine - Official Clinical Overview & Datasheet

### EXECUTIVE SUMMARY

The SHR (Super Hair Removal) machine represents a paradigm shift in phototrichology, positioning itself as the gold standard for high-volume medical aesthetics clinics and dermatology centers. Unlike traditional IPL or single-pulse diode lasers that rely on high fluence and thermal coagulation, the SHR modality leverages low-fluence, high-frequency cumulative heating to induce anagen follicle apoptosis without epidermal shock. This clinical architecture directly addresses the two largest barriers to adoption: patient pain perception and Fitzpatrick IV-VI treatment risk. The device's primary value proposition is painless, high-efficiency hair reduction with zero downtime, delivering superior return on investment through increased patient throughput and reduced consumable costs.



## CLINICAL ARCHITECTURE & DESIGN

The system is engineered around a military-grade diode laser stack (imported from either Jenoptik or DILAS) delivering a true 808nm wavelength, selected for its maximal melanin absorption to deep follicle depth ratio. The handpiece integrates a triple-spectrum capability via mechanical filter switching (755nm / 808nm / 1064nm), allowing clinicians to selectively target superficial, mid-dermal, or deep vascular chromophores. For thermal management, a four-stage cascade cooling system is employed: (1) dual TEC (Thermoelectric Cooler) modules achieve subzero contact plate temperatures, (2) a synthetic sapphire window provides 93% light transmittance while actively conducting heat away from the stratum corneum, (3) a closed-loop water circulation system with a dual-chamber pump dissipates heat from the laser stack, and (4) internal high-RPM fans exhaust residual thermal load. The console housing is a

corrosion-resistant ABS-alloy blend with sealed optical pathways to prevent dust ingress during extended clinical operation.

#### KEY INDICATIONS & CAPABILITIES

- SHR Pile-Up Mode: Proprietary pulse stacking algorithm delivering 1-10Hz frequency ramping. Low fluence (3-10 J/cm<sup>2</sup>) combined with rapid repetition raises target tissue temperature gradually to 45-50 ° C, achieving non-coagulative thermal injury. This eliminates the sharp "snap" sensation associated with standard lasers.

- In-Motion Technology: Integrated contact-based motion sensor disables emission if the handpiece remains stationary or lifts off the skin. This prevents geometric overlap burns and forces sweeping technique, reducing treatment time by 40% compared to stamping methods.

- Smart Clinical Database & UI: 10.4-inch capacitive touchscreen with pre-programmed protocols for 9 body zones (upper lip, axillae, bikini, legs, back, chest) across 6 Fitzpatrick skin types. The software adjusts fluence, pulse width, and delay automatically based on user-input patient parameters.

- Dual Cooling Modes: Epidermal protection selectable between "Freeze Mode" (sub-zero contact cooling for sensitive areas) and "Warm Mode" (neutral 10°C for large surface area treatments), preventing both thermal damage and cryogenic epidermal adhesion.

- Large Spot Size & Repetition Rate: 12mm x 12mm to 15mm x 15mm square

spot geometry eliminates diffraction patterns. Repetition rate up to 10Hz enables a full leg treatment in under 18 minutes.

## COMPLIANCE & STANDARDS

This device is manufactured under ISO 13485:2016 certified quality management systems for medical devices. It has received Medical CE (Class IIb) marking under the Medical Device Regulation (EU) 2017/745. For the United States market, the device holds FDA 510(k) clearance for permanent hair reduction (OTC over-the-counter not applicable, professional use only). Additional certifications include IEC 60825-1:2014 for laser safety, ROHS compliance for hazardous substance restrictions, and CB Scheme certification for global market entry. Each unit undergoes a 72-hour burn-in test and optical power calibration traceable to NIST standards.

## TECHNICAL SPECIFICATIONS

The following parameters are standardized for the professional SHR platform. All values measured at the handpiece output aperture under full operational load.

Parameter	Specification
-----------	---------------

Laser Type / Wavelength	808nm Diode Laser (Standard) / Triple wavelength 755/808/1064nm (Optional)
Fluence (Energy Density)	1 - 50 J/cm <sup>2</sup> (Adjustable, step 0.5 J/cm <sup>2</sup> )
Pulse Width	1 - 360 ms (Adaptive to frequency mode)
Repetition Rate (Frequency)	1 - 10 Hz (SHR mode: 1-10Hz automatic ramping)
Spot Size	15 mm x 15 mm (Standard) / 12 mm x 12 mm (Precision tip)
Cooling System	TEC (-15 ° C to +10 ° C) + Sapphire Contact + Closed-Loop Water + Wind
Screen Interface	10.4 inch TFT LCD capacitive touchscreen, 1024x768 resolution
Power Requirements	110-240 VAC, 50/60 Hz, 1800VA (Peak)

