

# Selective Photothermolysis Architecture Reference Document: Portable Carbon Peel Laser Unit

## SELECTIVE PHOTOTHERMOLYSIS ARCHITECTURE REFERENCE DOCUMENT: PORTABLE CARBON PEEL LASER UNIT

### 1. MECHANISM PHILOSOPHY

The Portable Carbon Peel Laser Unit operationalizes the principle of selective photothermolysis via a dual-phase chromophore amplification cascade. In Phase One, a topical carbon-based suspension (nanoparticle lotion) is applied to the epidermis, where it passively diffuses into sebaceous follicles, intercellular crevices, and comedonal openings. In Phase Two, the Q-switched or long-pulsed neodymium:yttrium-aluminum-garnet (Nd:YAG) 1064 nm or 532 nm emission interacts with the exogenous carbon chromophore. The resultant rapid thermoelastic expansion generates a micro-acoustic shockwave, mechanically lysing the carbon particles while transferring kinetic energy to adjacent sebum, keratinous debris, and cellular membranes. This mechanism spares surrounding water and melanin absorbers, enabling a non-ablative epidermal clearance profile.



## 2. ENERGY TRANSMISSION DYNAMICS

The unit deploys a compact articulating arm or lightweight flexible waveguide terminating in a precision-machined, ergonomic handpiece. Energy density (fluence) is digitally regulated from  $0.5 \text{ J/cm}^2$  to  $5.0 \text{ J/cm}^2$  in  $0.1 \text{ J/cm}^2$  increments, with a nominal spot size of 3 mm to 6 mm (circular or hexagonal). Pulse duration is fixed in the nanosecond regime (5-20 ns for Q-switched variants) or selectable picosecond extension (350-450 ps) depending on configuration. Repetition rate is adjustable from 1 Hz to 10 Hz, enabling controlled sweeping motion. The system incorporates a real-time fluence feedback loop via a pyroelectric sensor, maintaining output stability within  $\pm 5\%$  over 8-hour operational cycles.

## 3. CLINICAL ADVANTAGES

- Non-ablative clearance: No frank wounding of stratum corneum, reducing post-inflammatory hyperpigmentation risk in Fitzpatrick skin types III-V.
- Dual-action debridement: Simultaneous exfoliation of surface debris and sebaceous gland expression.
- Minimally invasive alternative to microdermabrasion and chemical peels with comparable efficacy for comedonal acne and photoaging.
- Rapid treatment sessions: 15-20 minutes for full-face application, enabling high patient throughput.
- No consumable carbon supply required beyond standard topical suspension; handpiece lens is self-cleaning with integrated air purge.

#### 4. CERTIFICATIONS MATRIX

Parameter	Specification
Laser Type / Wavelength	Nd:YAG / 1064 nm (primary), 532 nm (optional)
Pulse Duration	≤ 10 ns (Q-switched) or 380 ps (Pico option)
Spot Size	3 mm, 4 mm, 5 mm (interchangeable caps)
Max Fluence (1064 nm)	5.0 J/cm <sup>2</sup>

Repetition Rate	1-10 Hz (adjustable)
Cooling System	Forced air + TEC contact cooling (0-5°C)
Dimensions (WxDxH)	320 x 240 x 260 mm
Weight	8.5 kg (includes battery backup)
Power Input	100-240 VAC, 50/60 Hz, 200 VA
Aiming Beam	635 nm, <1 mW, Class 2

#### Regulatory / Safety | Compliance Status

FDA 21 CFR 1040.10 | Class II (with CDRH variance)

IEC 60825-1:2014 | Class 4 Laser Product (certified)

ISO 13485:2016 | Full Quality Management System

MDR (EU) 2017/745 | CE 2797 (Class IIb)

IEC 60601-2-22:2019 | Surgical laser equipment safety

RoHS 3 (EU) 2015/863 | Compliant, lead-free solder

#### 5. EXACT SPECIFICATIONS

- Laser Medium: Nd:YAG, diode-pumped, air-cooled
- Wavelengths: 1064 nm (primary carbon peel); optional 532 nm (superficial pigmentation)

- Output Mode: Q-switched (standard) or Pico-second (optional module)
- Max Fluence: 5.0 J/cm<sup>2</sup> (1064 nm) / 2.5 J/cm<sup>2</sup> (532 nm)
- Spot Size: 3 mm, 4 mm, 5 mm (interchangeable handpiece caps)
- Pulse Width: ≤ 10 ns (Q-switched) or 380 ps (Pico option)
- Repetition Rate: 1-10 Hz, adjustable via footswitch or touchscreen
- Aiming Beam: 635 nm, <1 mW, Class 2 laser diode
- Cooling: Forced air + TEC (thermoelectric) handpiece tip cooling (0°C to 5°C contact)
- Dimensions (base unit): 320 mm (W) x 240 mm (D) x 260 mm (H)
- Weight: 8.5 kg (including integrated battery backup, 15 min runtime)
- Input Power: 100-240 VAC, 50/60 Hz, 200 VA

## 6. OPERATIONAL ROI SCHEMATIC

Treatment volume estimation (single operator, 8-hour clinical day):

- Setup & patient intake: 5 min per patient
- Carbon suspension application: 5 min
- Laser pass (full-face, 3 mm spot, 50% overlap): 12 min (approx. 1200 pulses at 2 Hz)
- Post-treatment cooling & removal: 5 min
- Total cycle time: 27 min → 17 patients per 8-hour shift (assuming 15 min buffer for room turnover).

Per-treatment consumable cost: USD \$2.50 (carbon lotion, disposable eye shields, tip covers). Average clinic charge per carbon peel: USD \$200–\$350. Gross margin per session: \$197.50 – \$347.50. Payback period at 12 sessions/week (2.5 days): 4-6 months.

