

Selective Photothermolysis Architecture Reference Document: Q-Switched
Nanosecond Laser Workstation

SELECTIVE PHOTOTHERMOLYSIS ARCHITECTURE REFERENCE DOCUMENT:
Q-SWITCHED NANOSECOND LASER WORKSTATION

EXECUTIVE SUMMARY

This document provides a comprehensive technical and clinical overview of the Q-Switched Nanosecond Laser Workstation, a premium medical aesthetic device engineered for the precise clearance of dermal and epidermal pigmented lesions, professional tattoo removal, and skin rejuvenation. The platform leverages nanosecond-domain pulse generation to induce selective photothermolysis, ensuring maximal target chromophore destruction with minimized thermal diffusion to surrounding tissues. Designed for high-volume dermatology clinics and medical spas, this workstation integrates advanced solid-state laser resonators, an articulated articulated arm delivery system, and a closed-loop contact cooling engine.



CLINICAL ARCHITECTURE & DESIGN

The device architecture prioritizes clinical efficacy and operational longevity. The optical core utilizes high-stability Nd:YAG laser crystals producing fundamental emission at 1064nm, with frequency-doubling optics generating the 532nm wavelength. The Q-switching mechanism (Pockels cell or saturable absorber) compresses the pulse energy into a 5-20 nanosecond window, generating peak powers exceeding 50 MW. The energy is transmitted through a seven-articulation, counter-balanced articulated arm featuring dielectric-coated, high-damage-threshold mirrors to a lightweight ergonomic handpiece. A high-brightness Xenon flashlamp pumps the laser cavity, supported by a recirculating deionized water cooling loop and thermoelectric chillers to maintain optical alignment and flashlamp lifetime.

KEY INDICATIONS & CAPABILITIES

The Q-Switched Nanosecond Laser Workstation is indicated for:

- Professional Tattoo Removal (multi-color, amateur, traumatic)
- Benign Pigmented Lesions (lentigos, ephelides, café au lait, nevi of Ota/Ito)
- Recalcitrant Melasma (adjunctive, low-fluence protocol)
- Acne Scars & Pores (resurfacing via fractional handpiece, optional)
- Photoaged Skin (diffuse erythema, dyschromia)
- Nevus of Ota clearance

CONTRAINDICATIONS: Active local or systemic infection, history of keloid scarring, photosensitivity disorders, anticoagulation therapy, pregnancy, recent tanning, or Fitzpatrick Skin Types V-VI for 532nm (high epidermal melanin risk).

COMPLIANCE & STANDARDS

The workstation is manufactured under ISO 13485:2016 certified facilities and complies with:

- IEC 60825-1:2014 (Class 4 Laser Product safety)
- FDA 21 CFR 1040.10 & 1040.11 (with variance)
- Medical Device Directive 93/42/EEC (CE 0197)
- MDR (EU) 2017/745 Class IIb
- RoHS & REACH environmental directives

Built-in safety features include: aperture shutters, emergency stop, key switch,

ready indicator, audible countdown, footswitch interlock, and real-time energy monitoring.

TECHNICAL SPECIFICATIONS

Parameter	Specification
Laser Type / Wavelengths	Q-Switched Nd:YAG / 1064 nm & 532 nm (frequency-doubled)
Pulse Width (nominal)	5 – 20 nanoseconds (FWHM, user-selectable)
Peak Power (1064 nm)	> 50 MW (at max fluence, 6 mm spot)
Maximum Fluence (1064 nm)	12 J/cm ² (adjustable in 0.1 J/cm ² increments)
Maximum Fluence (532 nm)	6 J/cm ² (adjustable in 0.1 J/cm ² increments)
Spot Size Range	1 mm – 8 mm (continuous zoom or fixed handpiece tips)
Repetition Rate	1 – 10 Hz (single shot to continuous)
Energy Stability	≤ ± 5% over 8 hours (after warm-up)
Cooling System	Closed-loop: TEC + Sapphire contact plate (+2°C to +8°C) + recirculating deionized water + fans

Delivery System	7-joint articulated arm, counter-balanced, 160 cm reach
Display / UI	10.4-inch resistive touchscreen, preset protocols for lesion/tattoo type and Fitzpatrick scale
Dimensions (W x D x H)	42 cm x 68 cm x 98 cm (console only)
Weight (console)	68 kg (150 lbs)
Electrical Requirements	200-240 VAC, 50/60 Hz, 15A, single-phase
Laser Class	Class 4 (IEC 60825-1)
Operating Environment	15 °C – 30 °C, 30% – 75% relative humidity (non-condensing)
Warranty (Standard)	Laser head: 24 months / System: 12 months parts & labor

CLINICAL PROTOCOLS – STANDARD TREATMENT PARAMETERS

TATTOO REMOVAL (BLACK/DARK BLUE): 1064nm, 5-8 J/cm², 5-10 ns, 2-4 mm spot, 4-6 week intervals. Endpoint: immediate whitening (frosting).

TATTOO REMOVAL (RED/ORANGE/YELLOW): 532nm, 2-4 J/cm², 5-10 ns, 2-3 mm spot, 6-8 week intervals. Endpoint: frosting + mild purpura.

LENTIGINES (SOLAR LENTIGO): 532nm, 1.5-2.5 J/cm², 5-10 ns, 3-4 mm spot.

Endpoint: immediate darkening (no frosting required).

NEVUS OF OTA: 1064nm, 6-10 J/cm², 5-10 ns, 3-4 mm spot. Multiple sessions (3-8).

CLINICAL WORKFLOW & PATIENT MANAGEMENT

1. PRE-TREATMENT: Informed consent, skin type assessment (Fitzpatrick I-VI), test spot in discreet area (24-48h), topical anesthetic optional (30-60 min).
2. INTRA-TREATMENT: Don appropriate laser safety eyewear. Set fluence 10-15% below purpura threshold. Activate Sapphire cooling (preset to +4°C). Overlap <10%. Deliver pulses at 1-2 Hz.
3. POST-TREATMENT: Apply broad-spectrum SPF 50+ and occlusive dressing (if ablated). Instruct on cold compresses, avoid sun exposure, and no picking of micro-crusts. Adverse events: transient erythema, edema, purpura (532nm), hypo/hyperpigmentation (rare, usually resolves).

