

# Compressive Microvibration Device - Official Clinical Overview & Technical Datasheet

## EXECUTIVE SUMMARY

The Compressive Microvibration Device (CMD) represents a paradigm shift in non-invasive mechanotransduction therapy for aesthetic and rehabilitative medicine. Unlike thermal or photonic modalities, the CMD delivers precisely controlled, low-amplitude, high-frequency mechanical compression cycles to dermal and subdermal tissues. This proprietary actuation mechanism stimulates fibroblast activity, enhances local microcirculation, and initiates a cascade of regenerative extracellular matrix (ECM) remodeling without thermal injury or patient downtime. The system is engineered for deployment in medical spas, dermatology clinics, and physical medicine centers seeking an evidence-based, painless intervention for skin laxity, periorbital edema, and post-procedural recovery.



## CLINICAL ARCHITECTURE & DESIGN

The device platform integrates a closed-loop electromechanical actuator housed within a medical-grade, passively cooled chassis. Key design elements include:

- ACTUATOR TYPE: Balanced eccentric rotary mass oscillator with active frequency stabilization
- DELIVERY MODE: Unidirectional compressive wave train at predefined surface contact pressure
- FEEDBACK CONTROL: Piezoresistive force sensor array (8 points per cm<sup>2</sup>) for real-amplitude compliance
- OPERATOR INTERFACE: 10.1-inch capacitive touchscreen with haptic confirmation and treatment protocol library

- PATIENT SAFETY: Automatic shut-off on overtemperature (>42 ° C skin interface) or abnormal impedance

The system operates on a universal 100-240V AC input, converting to a 48V DC actuator drive. The handpiece is autoclave-capable for multi-patient clinical environments.

#### KEY INDICATIONS & CAPABILITIES

- PRIMARY INDICATIONS: Mild-to-moderate facial and neck skin laxity, orbital puffiness, post-liposuction fibrosis modulation, and acceleration of ecchymosis resolution.
- SECONDARY APPLICATIONS: Pre-treatment tissue preconditioning (enhances topical penetration) and post-laser edema reduction.
- CONTRAINDICATIONS: Active malignancy, metallic implants in treatment field, pregnancy, acute thrombophlebitis.
- TREATMENT PROFILE: Painless (mean VAS 0.4/10), no anesthetic required, zero thermal injury, no purpura.

#### COMPLIANCE & STANDARDS

The Compressive Microvibration Device is manufactured under ISO 13485:2016

certified quality management systems. It complies with:

- IEC 60601-1 (Medical electrical equipment – general safety)
- IEC 60601-1-2 (Electromagnetic compatibility – emissions and immunity)
- Classification: Class IIa (Medical Device Regulation (EU) 2017/745)
- FDA 510(k) cleared for adjunctive treatment of superficial soft tissue (Kxxxxxx  
– reference available on request)

<b>Parameter</b>	<b>Specification</b>
Actuation Type	Electromechanical eccentric rotary mass oscillator
Frequency Range	20–60 Hz ( $\pm 2\%$ stability)
Peak Amplitude	0.2 – 1.2 mm (programmable step 0.05 mm)
Contact Pressure Sensing	Piezoresistive array, 0–10 N range
Power Supply	100–240 V AC, 50/60 Hz, 120 VA max
Dimensions (Base Unit)	280 x 260 x 180 mm (W x D x H)
Weight (Base Unit)	4.2 kg
Handpiece Mass	280 g (without cable)
Handpiece Cable Length	2.5 m (shielded)
Standby Noise	< 35 dB(A)
Operating Noise	< 55 dB(A) at 50 Hz

Display	10.1" TFT LCD, 1280x800, capacitive touch
Treatment Memory	500 patient records, exportable via USB-C
Certifications	CE MDR (Class IIa), FDA 510(k), ISO 13485

## CLINICAL PROTOCOLS

### STANDARD PROTOCOL (FACIAL SKIN LAXITY)

- Frequency: 30–50 Hz (titratable)
- Amplitude: 0.5–1.2 mm (peak-to-peak displacement)
- Contact force: 2–5 N (feedback-controlled)
- Session duration: 10–15 minutes per zone
- Recommended course: 6–8 sessions, 1x weekly

### POST-LASER RECOVERY PROTOCOL

- Frequency: 20–30 Hz (low mechanical stress)
- Amplitude: 0.3–0.6 mm
- Session duration: 5–8 minutes
- Timing: Immediately post-ablative treatment and 24 hours later

## PERIORBITAL EDEMA PROTOCOL

- Frequency: 40–60 Hz
- Amplitude: 0.2–0.4 mm (reduced due to thin orbital skin)
- Contact force: 1–2 N
- Session duration: 5 minutes per orbit

All protocols are pre-programmed into the device library, with manual override for advanced clinicians.

