

Liftera Alternative Line HIFU Device - Clinical Architecture & Performance
Reference Manual

CLINICAL ARCHITECTURE & PERFORMANCE REFERENCE MANUAL: LIFTERA
ALTERNATIVE LINE HIFU DEVICE

EXECUTIVE SUMMARY

The Liftera Alternative Line HIFU Device represents a paradigm shift in non-invasive skin laxity treatment, integrating advanced High-Intensity Focused Ultrasound (HIFU) technology with an ergonomic, high-efficiency platform. Engineered for the modern dermatology clinic and medical spa, this system delivers precision micro-thermal coagulation zones at predetermined depths, stimulating neocollagenesis and elastin regeneration without compromising the epidermal layer. This document provides a comprehensive technical overview, detailing the clinical architecture, performance specifications, and operational protocols essential for integrating this premium asset into a high-ROI aesthetic practice. The device is distinguished by its exceptional energy delivery consistency, advanced cooling mechanisms for patient comfort, and an intuitive user interface designed for rapid treatment workflows.



CLINICAL ARCHITECTURE & DESIGN

The foundational architecture of the Liftera Alternative Line is built upon a robust, high-frequency ultrasound generator that converts electrical energy into precise acoustic energy. This energy is focused through a proprietary transducer housed within the treatment handpiece, creating a converged beam that targets the superficial musculoaponeurotic system (SMAS) and deep dermal layers. The system employs a sophisticated beam-forming algorithm to ensure the focal point remains constant, delivering consistent energy fluence regardless of operator hand speed. The chassis is constructed from a durable, medical-grade polymer alloy, providing thermal management and structural integrity for high-volume clinical use. Internal subsystems include a redundant temperature monitoring network and a high-efficiency switching power supply that maintains stable output across varying line voltages, ensuring global

operational reliability.

KEY INDICATIONS & CAPABILITIES

The Liftera Alternative Line is indicated for the non-invasive lifting and tightening of facial and submental skin. Its primary capabilities include the reduction of mild to moderate skin laxity, the improvement of jawline definition, and the elevation of brow position. The device is uniquely equipped to treat a broad spectrum of phototypes due to its non-ablative, non-thermal epidermal mechanism. It is capable of creating targeted thermal lesions at depths of 1.5mm, 3.0mm, and 4.5mm, allowing for customized treatment planning based on individual patient anatomy and skin thickness. This multi-depth capability ensures comprehensive dermal remodeling, from the superficial papillary dermis to the deep reticular dermis and SMAS layer, providing a versatile solution for facial rejuvenation and contouring.

COMPLIANCE & STANDARDS

The Liftera Alternative Line HIFU Device is manufactured in adherence to the highest international quality and safety standards. The system conforms to the essential requirements of the Medical Device Directive (MDD) 93/42/EEC and is CE marked for distribution within the European Economic Area. Additionally,

the manufacturing facility is certified under ISO 13485:2016 for quality management systems specific to medical devices. The device has undergone rigorous electromagnetic compatibility (EMC) testing in accordance with IEC 60601-1-2 and electrical safety testing per IEC 60601-1. These certifications guarantee that the equipment is safe for both the patient and the operator when used in a clinical environment. The unit also includes integrated safety interlocks that prevent inadvertent energy discharge, ensuring compliance with global laser and ultrasound safety standards.

TECHNICAL SPECIFICATIONS

The operational performance of the Liftera Alternative Line is defined by a precise set of metrics designed to ensure predictable and repeatable clinical outcomes. The system parameters are configurable via a high-resolution touchscreen interface, allowing clinicians to select the optimal energy and depth settings for each unique patient presentation. The table below encapsulates the core performance parameters and physical specifications of the device, providing a definitive reference for clinical engineers and procurement specialists.

Parameter	Specification
Device Type	High-Intensity Focused Ultrasound

	(HIFU)
Focal Depths	1.5mm, 3.0mm, 4.5mm (selectable)
Energy per Shot	0.05 J to 0.35 J (adjustable in 0.01 J increments)
Number of Transducers per Cartridge	7
Cartridge Shot Count	500 shots per cartridge
Treatment Spacing	1.0mm to 3.0mm (user-selectable)
Display Interface	10.1-inch multi-touch capacitive touchscreen
Cooling System	Advanced forced-air and integrated thermoelectric (TEC) cooling
Dimensions (Main Unit)	380mm (W) x 450mm (D) x 320mm (H)
Weight (Main Unit)	12.5 kg
Power Supply	100-240V AC, 50/60Hz, 2.5A
Operational Environment	Temperature: 10°C to 30°C; Humidity: 30% to 75% RH

CLINICAL PROTOCOLS

Optimized clinical application of the Liftera Alternative Line requires adherence to a structured treatment protocol to maximize efficacy and ensure patient

safety. The recommended protocol begins with a thorough patient consultation and skin assessment to identify suitable treatment zones. Pre-treatment, the targeted area must be cleansed and degreased, and a coupling gel applied to ensure optimal acoustic transmission. The handpiece is positioned perpendicular to the skin surface, and energy is delivered in a grid-like pattern, with precise spacing between shots to create a contiguous treatment field. The integrated, dynamic cooling system automatically manages surface temperature, ensuring patient comfort throughout the procedure. The post-treatment regimen includes the application of a hydrating serum and sunblock to support the skin's natural regenerative process, with patients typically experiencing mild erythema that resolves within a few hours.

