

PDF Document

```
```json
```

```
{
```

```
  "title": "Cryo Thermal Contrast Therapy System - Official Clinical Overview &  
Technical Datasheet",
```

```
  "content": "CRYO THERMAL CONTRAST THERAPY SYSTEM\nOFFICIAL  
CLINICAL OVERVIEW & TECHNICAL DATASHEET\n\nEXECUTIVE  
SUMMARY\n\nThis document provides a comprehensive technical and clinical  
overview of the Cryo Thermal Contrast Therapy System, a premium,  
multi-modality medical aesthetic device. Designed for deployment in high-end  
med spas, dermatology clinics, and rehabilitation facilities, this system  
integrates advanced cryotherapy, thermotherapy, and rapid contrast therapy  
(RCT) within a single, intelligent platform . The device leverages a proprietary  
dual-reservoir thermal engine and a smart user interface to deliver precise,  
controlled thermal energy, enhancing patient comfort, accelerating recovery,  
and expanding clinical treatment capabilities. This datasheet outlines the  
system's clinical architecture, technical specifications, compliance standards,  
and operational protocols, establishing its role as a versatile asset in modern  
aesthetic and therapeutic medicine.\n\n[IMAGE_1]\n\nCLINICAL  
ARCHITECTURE & DESIGN\n\nThe Cryo Thermal Contrast Therapy System is  
engineered with a sophisticated internal architecture centered on a  
dual-reservoir heating and cooling system. This design overcomes the
```

limitations of single-reservoir devices by providing dedicated, pre-conditioned thermal sources, allowing for nearly instantaneous switching between hot and cold modalities. Each reservoir is served by an independent energy source and flow circuit, complete with pressure relief and check valves to ensure operational integrity and precise temperature control. A high-performance Thermoelectric Cooler (TEC) module serves as the primary thermal engine, offering accurate temperature modulation without the need for ice or external refrigerants. Electronically controlled three-way valves ensure that only the selected thermal medium circulates to the therapy handpiece at any given time, enabling seamless transitions during contrast therapy cycles. The system incorporates a closed-loop biofeedback mechanism via embedded skin temperature sensors, continuously monitoring the treatment site to maintain the desired therapeutic temperature and enhance patient safety.

**KEY INDICATIONS & CAPABILITIES**

This system is indicated for a broad spectrum of clinical applications, with the primary goal of reducing pain, inflammation, and edema associated with post-surgical and traumatic injuries. In the aesthetic context, it serves as a critical adjunct to energy-based procedures, effectively minimizing patient discomfort, erythema, and downtime. The system offers three primary treatment modalities:

1. **Cryotherapy (Cold Therapy):** Induces vasoconstriction to reduce blood flow, inflammation, and muscle spasms. It effectively numbs the treatment area, providing a potent analgesic effect that is crucial for patient comfort during and after invasive

aesthetic procedures.\n2. **Thermotherapy (Heat Therapy):** Promotes vasodilation to increase circulation, support the lymphatic system, and stimulate the body's natural healing processes. It is ideal for managing chronic pain, stiffness, and preparing tissues for treatment.\n3. **Contrast Therapy (Alternating Hot and Cold):** Rapidly cycles between hot and cold to create a 'pumping' action in the vasculature. This process increases blood flow to the injury site without causing additional swelling, significantly accelerating the removal of metabolic waste and the delivery of oxygen and nutrients, thereby reducing recovery time .\n\nCOMPLIANCE & STANDARDS\n\nThe Cryo Thermal Contrast Therapy System is manufactured in adherence to the highest standards of medical device quality and safety. It is FDA-cleared as a Class II medical device and complies with all applicable standards for electromagnetic compatibility (EMC) and electrical safety, including IEC 60601 series. The system incorporates multiple built-in safety features, including over-temperature protection, automatic shut-off timers, and continuous pressure monitoring to prevent adverse events . The software interface is validated to ensure accurate program execution and data integrity, making it a reliable and compliant solution for clinical environments. The system is intended for use by or on the order of licensed healthcare professionals in rehabilitation facilities, outpatient clinics, and athletic training settings .\n\nTECHNICAL SPECIFICATIONS\n\n[TABLE\_1]\n\nCLINICAL PROTOCOLS\n\nEffective use of the Cryo Thermal Contrast Therapy System requires adherence to established

clinical protocols tailored to the specific treatment and patient condition. The device features an intuitive touchscreen interface that provides access to pre-set, evidence-based treatment options for pain management, post-operative recovery, and aesthetic procedure support . Clinicians can also utilize the system's advanced software to create and save customized protocols, adjusting parameters such as temperature setpoints, treatment duration, and contrast cycle timing to meet individual patient needs. For contrast therapy, typical cycles alternate between 5 minutes of cold and 5 minutes of heat over a 20-minute session . During cryotherapy, the system's integrated biofeedback loop maintains precise skin contact temperature, which is critical for ensuring patient safety and achieving optimal therapeutic outcomes. It is recommended that the treatment area be clean and dry before application, and the therapy wrap secured to ensure uniform contact .\n\n[IMAGE\_2]",

"images": [

"A high-quality 4K realistic promotional image showing a sleek, modern generic aesthetic machine in a bright, luxurious clinical setting, high tech medical vibe, unbranded chassis, no text.",

"A high-quality 4K realistic close-up image showing the premium generic treatment handpiece, sapphire cooling tip, or smart interface of the aesthetic equipment, professional studio lighting, unbranded components, no text."

],

"tables": [

```
{  
  "headers": ["Parameter", "Specification"],  
  "rows": [  
    ["Cooling Range", "34°F - 59°F (1°C - 15°C) "],  
    ["Heating Range", "Up to 113°F (45°C) "],  
    ["Reservoir System", "Dual Temperature, Dual Reservoir "],  
    ["Cooling Technology", "Thermoelectric Cooler (TEC) / Peltier Cell "],  
    ["Max Treatment Time (Cold)", "90 Minutes "],  
    ["Max Treatment Time (Heat)", "30 Minutes "],  
    ["Interface", "Touchscreen Computer Interface "],  
    ["Safety Compliance", "FDA Class II Cleared , IEC 60601"],  
    ["Wraps Supported", "Up to 4 Thermal Wraps (Model Dependent) "  
  ]  
}  
]  
}
```